



**MILLENNIUM  
SURGICAL**




SURGICALINSTRUMENTS.COM

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Hand Tables & Retractors

## Hand Tables & Retractor Components



18-TH110 and 18-TH111 are not included in the child and infant sets. They are sold separately.

	Quantity	Cat. No.	Description
<b>A</b>	1 ADULT	18-TH101	The <b>tables</b> maintain enough weight to resist all passive torsion force of the pronated or supinated forearm. Around its perimeter numerous notches provide anchor spots for the ball chains. Slots serve as placements for the rubber bands' sliding hooks. The tables can be used facing either way to accommodate both left and right hand surgical procedures.
<b>B</b>	1 CHILD	18-TH201	
<b>C</b>	1 INFANT	18-TH301	
<b>D</b>	1 set of 2 pcs.	18-TH102	<b>Elevated chain holders</b> are positioned anywhere around the edge of the table and serve as elevated anchor points for chains wherever an elevated pull or clearance above the flexed fingers is needed.
<b>E</b>	1 set of 2 pcs.	18-TH103	<b>Chain holder hooks</b> enable one chain to provide an anchor point for another chain. They are used to facilitate proximal retraction along the axis of the limb (See picture).
<b>F</b>	1 set of 2 pcs.	18-TH104	<b>Ball chains with hooks</b> , 32cm long. <b>Ball chains with hooks</b> , 21cm long. Both chains are used for wound edge retraction and retraction of fingertips. The chain, engaged in the table edge, or the elevated chain holders, or in the chain holder hooks, has an interval of 6cm.
<b>G</b>	1 set of 4 pcs.	18-TH105	
<b>H</b>	1 set of 2 pcs.	18-TH106	<b>Ball chains without hooks</b> , 26cm long, are used to secure arm and wrist fixation. They are secured at each end to the table edge.
<b>I</b>	1 set of 4 pcs.	18-TH107	<b>Tendon hooks</b> , when attached to the hooks of TH-104/105, facilitate the traumatic retraction of tendons and wound edges.
<b>J</b>	1 set of 4 pcs.	18-TH108	<b>Rubber band sliding hooks</b> are placed at any point along the slots to provide anchor points for the rubber bands.
<b>K</b>	1 set of 5 pcs.	18-TH109	<b>Rubber bands</b> are used to provide a downward pull to the various digits at any point along their length.
<b>L</b>	1	18-TH110	<b>Stainless Steel Carrying Case</b>
<b>M</b>	1	18-TH111	<b>Oval box</b> 130 x 8mm, for parts: TH-102 through TH-109.



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Traction Tower



100FT Traction Device

Autoclavable stainless steel construction with a telescoping traction post. This quick-mounting unit, combined with the finger traps, is the perfect system for traction procedures (weights sold separately).



EQ350 Equalizing Device

Provides equal traction to as many as five digits and is autoclavable.



TT200 Traction Tower

This simple but effective unit mounts with ease to the O.R. table accessory rail with a clamp which is attached to the vertical post. This 36" rod is height-adjustable using the clamp. The 30" boom angle can be adjusted and rotated 360°. This length provides the versatility to place the tower on either side of the O.R. table.



**MILLENNIUM  
SURGICAL**

**SURGICALINSTRUMENTS.COM**

UNMATCHED SERVICE · UNLIMITED INSTRUMENTS

822 Montgomery Ave #205  
Narberth, PA 19072  
(800) 600-0428

# Traction Tower



## Stainless Steel Finger Traps

These autoclavable single wire traps can be ordered in (4) sizes.

Single Mode  
100XSW Pediatric  
200SW Small  
300MW Medium  
400LW Large



## Autoclavable Finger Traps

Fabric mesh autoclavable finger traps in five sizes. Available in single and double modes.

	Single Mode	Double Mode
Pediatric	100XS	102XS
Small	200S	202S
Medium	300M	302M
Large	400L	402L
X-Large	500XL	502XL